

Asepti Neutral Wipes

Neutral Detergent Wipes - 325 x Pail, 100 x Soft Pack

ARTG 239652



Features

- Asepti Neutral Wipes can be used on most hard surfaces, including to clean medical devices and other therapeutic goods, effectively removing common soils. They are suitable for use in hospitals, doctors and dental surgeries and other healthcare facilities.
- Each Asepti Neutral Wipe cloth is moist and of a generous size providing good coverage. The fabric is a thermally bonded viscose ensuring a soft, weighty and absorbent wipe perfect for soil removal.
- Asepti Neutral Wipes are available in a variety of handy pack sizes: large pail for heavy use areas (Operating Theatres, CSSDs etc.) and soft packs for lower use areas.
- Asepti Neutral Wipes are disinfectant free and are conveniently packed in 100 wipe flat packs and 325 wipe pails.
- Brackets are available to hold the product in place (on the wall, trolley etc.) which enables the product to be located within easy reach.

Product Code	Product Description	Pack Size	Quote
7100750	Asepti Neutral Wipes	3 x 325 wipes, Pail	
7100848	Asepti Neutral Wipes	12 x 100 wipes, Flat Pk	
9996-1208, 10112033	Wall Bracket to suit Pail	ea	
99001137, 10111807	Plastic Holder to suit Flat Pack	ea	
99960090, 10111973	Wire Holder to suit Flat Pack	ea	

Safety Information

Do not mix with any other products. For further information and safety instructions see product label and Safety Data

Sheet. ALWAYS READ THE LABEL, FOLLOW DIRECTIONS FOR USE.

Directions For Use

- For best results apply to dry surface in a normal wiping motion.
- When cleaning blood or body soils from hard surfaces, wear protective gloves (not required when wiping surfaces which are free of visible soils).
- Dispose of wipe and glove after each use and wash hands with the appropriate [Ecolab Endure hand wash](#).
- Safe for use on most surfaces. Can be used to clean medical devices and other therapeutic goods.
- Not suitable for use as a baby wipe. Wipes are single use only.

Related Products

- [Asepti Active Wipes](#)
- [Asepti Alcohol Wipes](#)
- [Endure Hand Wash](#)

Applications

- Theatre
- Dental
- Laboratory
- General Medical
- Patient Room
- Emergency
- Waiting Rooms
- Office / Meeting Rooms
- Endoscopy
- Residential Aged Care
- Life Science
- Healthcare Facilities
- Wards
- Day Surgeries
- Veterinary

Resources

- Safety Data Sheet (SDS)
- Brochure