


Guardian Clean

Washer disinfectant tablets for laboratory utensils and glassware



Features

- Guardian Clean tablets are designed for automated cleaning of laboratory utensils and glassware in food science, environmental laboratories and hospital pantry washers.
- Guardian Clean tablet is a dual enzymatic, non-caustic, non-chlorinated, non-phosphate formula with biodegradable surfactants to remove even the toughest soils in one pass.
- Guardian Clean has been shown to assist with reducing commonly occurring pathogens Escherichia coli.

Product Code	Product Description	Pack Size	Quote
9046400	Guardian Clean Tablets, 4 Kg / 200 Tablets	ea	

Characteristics

- Appearance: Tablet
- Colour: Opaque, white with coloured particles
- pH: 9.0 - 11.0, (1%)
- Odour: Slight

Safety Information

Do not mix with any other products. Further information, safety instructions and health warnings, see product label and Safety Data Sheet.
 ALWAYS READ THE LABEL. FOLLOW DIRECTIONS OF USE.

Directions For Use

- Check product label for detailed use information.
- Take one tablet of Guardian Clean Tablet and open the plastic foil.
- Place the tablet in the product reservoir and start normal program.
- Use at 50°C or lower for optimum cleaning results. Use concentration = 1 tab per 8.5 litres of water.

For professional use only

- This detergent is ideally suited to under bench machines. Select low temperature washing cycles, wash full loads and do not exceed the recommended dosage. This will minimise both energy and water consumption without compromising sanitation result.
- Recommended storage temperature between 0°C - 40°C in original container only.

Relevant Product

- [Smeg GW0260](#)
- [Smeg GW2145](#)
- [Smeg GW4060](#)
- [Smeg GW4190](#)
- [Smeg GW3060](#)
- [Smeg GW4090](#)
- [Smeg GW6290](#)
- [Smeg GW6010](#)
- [Smeg GW7010](#)
- [Smeg GW7015](#)
- [Belimed WD290](#)

Applications

- Life Science
- Laboratory Infection Control
- Healthcare Facilities
- Central Sterile Departments

Resources

- safety data sheet (sds)
- Video
- Brochure