


Smartline SmartSink











Electric Height Adjustable Washing Stations improve productivity





Features

- Smartline SmartSink offers your reprocessing department an electric height adjustable washing station which improves productivity and reduces the risk of staff injury.
- Smartline SmartSinks are manufactured in Australia and constructed out of 304 grade stainless steel.
- The Smartline range of SmartSinks offer a complete, customisable washing station solution for your departments requirements.
- Smartline SmartSinks feature a 300mm high tapered splashback and tapered benchtop with wetdeck perimeter, offering smooth surfaces for easy cleaning for infection prevention.
- SmartSinks include easy-to-use electric height adjustment with anti-collision technology from 850mm to 1080mm off the ground. Controls are positioned under the sink top offering the user full control and easy access.
- Smartline SmartSinks offer so many customisable options; multiple bench surface, two different or customisable bowl depth, rectangular, circular, single and/or double sink options, sink inserts, shelves, trays and racks and tapware to name a few.
- All Smartline SmartSinks are manufactured with adjustable bolt down feet.

Product Code	Product Description	Pack Size	Quote
SMT-11.SS-135-SR, 10118352	Smartline 1350mm SmartSink Single Bowl - Rectangular	ea	
SMT-11.SS-175-SR, 10118356	Smartline 1750mm SmartSink Single Bowl - Rectangular	ea	

Product Code	Product Description	Pack Size	Quote
SMT-11.SS-195-SR, 10118360	Smartline 1950mm SmartSink Single Bowl - Rectangular	ea	
SMT-11.SS-220-SR, 10118365	Smartline 2200mm SmartSink Single Bowl - Rectangular	ea	
SMT-11.SS-260-SR, 10118368	Smartline 2600mm SmartSink Single Bowl - Rectangular	ea	
SMT-11.SS-320-SR, 10118372	Smartline 3200mm SmartSink Single Bowl - Rectangular	ea	
SMT-11.SS-175-DR-1D, 10118354	Smartline 1750mm SmartSink Double Bowl - Rectangular	ea	
SMT-11.SS-195-DR-1D, 10118358	Smartline 1950mm SmartSink Double Bowl - Rectangular	ea	
SMT-11.SS-220-DR-1D, 10118362	Smartline 2200mm SmartSink Double Bowl - Rectangular	ea	
SMT-11.SS-220-DR-2D, 10118363	Smartline 2200mm SmartSink Double Bowl - Rectangular	ea	
SMT-11.SS-260-DR-1D, 10118367	Smartline 2600mm SmartSink Double Bowl - Rectangular	ea	
SMT-11.SS-320-DR-1D, 10118371	Smartline 3200mm SmartSink Double Bowl - Rectangular	ea	
SMT-11.SS-360-DR-1D, 10118374	Smartline 3600mm SmartSink Double Bowl - Rectangular	ea	
SMT-11.SS-UT13560, 10118381	Smartline Under Sink Tray - suit 1350x600 Sink Units	ea	
SMT-11.SS-UT13570, 10118382	Smartline Under Sink Tray - suit 1350x700 Sink Units	ea	
SMT-11.SS-UT17560, 10118383	Smartline Under Sink Tray - suit 1750x600 Sink Units	ea	
SMT-11.SS-UT17570, 10118384	Smartline Under Sink Tray - suit 1750x700 Sink Units	ea	
SMT-11.SS-UT19570, 10118386	Smartline Under Sink Tray - suit 1950x700 Sink Units	ea	
SMT-11.SS-UT22060, 10118387	Smartline Under Sink Tray - suit 2200x600 Sink Units	ea	
SMT-11.SS-UT22070, 10118388	Smartline Under Sink Tray - suit 2200x700 Sink Units	ea	
SMT-11.SBT2080, 10118349	Smartline SmartSink Splashback Tray 20x80cm	ea	
SMT-11.SS-UT26060, 10118389	Smartline Under Sink Tray - suit 2600x600 Sink Units	ea	
SMT-11.SS-UT26070, 10118390	Smartline Under Sink Tray - suit 2600x700 Sink Units	ea	

Product Code	Product Description	Pack Size	Quote
SMT-11.SBS2050, 10118346	Smartline Smartsink Splash Back Shelf - 20x50cm	ea	
SMT-11.SBS2080, 10118347	Smartline Smartsink Splash Back Shelf - 20x80cm	ea	
SMT-11.SBT2050, 10118348	Smartline SmartSink Splashback Tray 20x50cm	ea	
SMT-11.SBT2080, 10118349	Smartline SmartSink Splashback Tray 20x80cm	ea	
SMT-11.SS-CHEMSHLF, 10118377	Smartline SmartSink Chemical Bottle Shelf	ea	
SMT-11.SBCB, 10118345	Smartline SmartSink Clamp Mount for Accessories	ea	
SMT-51.TRYRCK-1X5, 10118577	Smartline SmartSink Stainless Steel Tray Rack for Autoplas Trays	ea	
SMT-51.TRDIN-1X5, 10118576	Smartline SmartSink Stainless Steel Tray Rack for DIN Trays	ea	
SMT-11.SS-GF152510-SS-SH, 10118378	Smartline SmartSink Pre-Rinse Unit, 10 inch spout	ea	
SMT-11.TAP-CL10010, 10118393	Smartline SmartSink Lever Mixing Set	ea	
SMT-11.TAP-CL15035, 10118394	Smartline SmartSink Lever Sink Set	ea	
SMT-11.TAP-JETF3000W, 10118397	Smartline SmartSink Pre-Rinse Unit	ea	
SMT-11.TAP-DURA295, 10118395	Smartline Tapware - Laboratory Stop Tap - Hob Mounted Recirculating	ea	
SMT-11.TAP-DURA299, 10118396	Smartline Tapware - Laboratory Stop Tap - Hob Mounted	ea	
SMT-47.DINA-JET, 10118512	Smartline Dina Cleaning Jet Pistol with Standard Set	ea	
SMT-47.HELTA-JET, 10118513	Smartline Hela Cleaning Jet Pistol with Standard Set, Stainless Steel	ea	

Relevant Product

- [Smartline Benches](#)
- [Aniosyme Foam](#)
- [Asepti Clean](#)
- [Asepti Epizyme Ultra](#)
- [Asepti Finish](#)
- [Asepti Metal Bright](#)
- [Asepti Release Plus](#)
- [Asepti Wash](#)
- [OptiPro Multi Enzyme](#)
- [Optipro Pre-Treatment Gel](#)
- [Sekusept Multi-Enzyme](#)
- [Water Testing, Table 7.3](#)

- [Asepti Multizyme](#)
- [Asepti Neutrazyme](#)
- [Asepti Newgen](#)
- [Asepti Rapidzyme PR](#)
- [Central Service Program](#)
- [Transportation Trolleys](#)

Applications

- Endoscopy
- Central Sterile Department
- Healthcare Facilities

Resources

- video
- Brochure